


Product Name: Recombinant Human FARP2 Protein
Catalog #: PEH8522



Summary

Name	FARP2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O94887
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	FRYSGRTQKQ LVDYFKDSGM KRIPYERRHS KTHTSVRALT ADLPKQSISF PEGLRTPASP SSANAFYSLs PSTLVPSGLP EFKDSSSSLT DPQVSYVKSP AAERRSGAVA GGPDTPSAQP LGPPALQPGP GLSTKSPQPS PSSRKSPSL SPAFQVPLGP AEQGSSPLLS PVLSDAGGAG

image



Background

Alternative Names	FERM, ARHGEF and pleckstrin domain-containing protein 2 (FERM domain-including RhoGEF) (FIR) (FERM, RhoGEF and pleckstrin domain-containing protein 2) (Pleckstrin homology domain-containing family C member 3) (PH domain-containing family C member 3)
--------------------------	---

Product Name: Recombinant Human FARP2 Protein
Catalog #: PEH8522



Note

For Research Use Only , Not for Diagnostic Use