


Product Name: Recombinant Human FOXO6 Protein
Catalog #: PEH8524



Summary

Name	FOXO6_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	A8MYZ6
Host	E.coli
Species	Human
Molecular Mass	51 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MAAKLRAHQV DVDPDFAPQS RPRSCTWPLP QPDLAGDEDG ALGAGVAEGA EDCGPERRAT APAMAPAPPL GAEVGPLRKA KSSRRNAWGN LSYADLITKA IESAPDKRLT LSQIYDWMVR YVPYFKDKGD SNSSAGWKNS IRHNLSLHTR FIRVQNEG TG KSSWWMLNPE GGKTGKTPRR RAVSMDNGAK FLRIKGKASK KKQLQAPERS PDDSSPSAPA PGPVPAAAKW AASPASHASD DYEAWADFRG GGRPLLGEAA ELEDDEALEA LAPSSPLMYP SPASALSPAL GSRCPGELPR LAELGGPLGL HGGGGAGLPE GLLDGAQDAY GPRPAPRPGP VLGAPGELAL AGAAAAYPGK GAAPYAPPAP SRSALAHPI S LMTLPGEAGA AGLAPPGHAA AFGGPPGGLL LDALPGPYAA AAAGPLGAAP DRFPADLDLD MFGSLECDV ESIILNDFMD SDEMDFNFDS ALPPPPPGLA GAPPPNQSWV PG

image



Background

Product Name: Recombinant Human FOXO6 Protein
Catalog #: PEH8524



Alternative Names

Forkhead box protein O6

Note

For Research Use Only , Not for Diagnostic Use