



Product Name: Recombinant Human FMN1 Protein
Catalog #: PEH8525



Summary

Name	FMN1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q68DA7
Host	E.coli
Species	Human
Molecular Mass	43 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SLSQEPDEHP GDIFFKQTPT KDILTELYKL TTERERLLTN LLSSDHILGI TMGNQEGKLQ ELSVSLAPED DCFQSAGDWQ GELPVGPLNK RSTHGNKKPR RSSGRRESFG ALPQKRTRRK GRGGRESAPL MGKDKICSSH SLPLSRTRPN LWVLEEKGNL LPNGALACSL QRRESCPPDI PKTPDIDLGF GSFETAFKDT GLGREVLPPD CSSTEAGGDG IRRPPSGLEH QQTGLSESHQ DPEKHPEAEK DEMEKPAKRT CKQKPVSKVV AKVQDLSSQV QRVVKTHSKG KETIAIRPAA HAEFVPKADL LTLPGAEGA HGSRRQGKER QGDRSSQSPA GETASISSVS ASAEGAVNKV PLKVIESEKL DEAPEGKRLG

image



Background

Alternative Names Formin-1 (Limb deformity protein homolog)

Product Name: Recombinant Human FMN1 Protein
Catalog #: PEH8525



Note

For Research Use Only , Not for Diagnostic Use