


Product Name: Recombinant Human FGD1 Protein
Catalog #: PEH8526



Summary

| | |
|------------------------------|---|
| Name | FGD1_HUMAN |
| Purity | Greater than 85% as determined by reducing SDS-PAGE |
| Endotoxin level | Please contact with the lab for this information |
| Tag | C-His |
| Accession # | P98174 |
| Host | E.coli |
| Species | Human |
| Molecular Mass | 13 kDa |
| Buffer | Phosphate buffered saline |
| Shipping | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below. |
| Stability&Storage | Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| Sequence | PSPQHHRALR FSYHLEGSQP RPLHQGNRI LVKSLSLDPG QSLEPHPEGP QRLRSDPGPP TETPSQRPS LKRAPGPKPQ VPPKPSYLQM PRMPPPLEPI |

image



Background

Alternative Names FYVE, RhoGEF and PH domain-containing protein 1 (Faciogenital dysplasia 1 protein) (Rho/Rac guanine nucleotide exchange factor FGD1) (Rho/Rac GEF) (Zinc finger FYVE domain-containing protein 3)

Note

For Research Use Only , Not for Diagnostic Use