



Product Name: Recombinant Human GPS2 Protein
Catalog #: PEH8527



Summary

Name	GPS2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q13227
Host	E.coli
Species	Human
Molecular Mass	19 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QEEEEVDKMM EQKMKEEQER RKKKEMEERM SLEETKEQIL KLEEKLLALQ EEKHQLFLQL KVLHHEEKR RRKEQSDLTT LTAAYQQSL TVHTGTHLLS MQGSPGGHNR PGLMAADRA KQMFQPQLT TRHYVGSAAA FAGTPEHGQF

image



Background

Alternative Names G protein pathway suppressor 2 (GPS-2)

Note

For Research Use Only , Not for Diagnostic Use