


Product Name: Recombinant Human GRSF1 Protein
Catalog #: PEH8528



Summary

Name	GRSF1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q12849
Host	E.coli
Species	Human
Molecular Mass	42 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PPPEYELAPS KLEEEVDDVF LIRAQGLPWS CTMEDVLNFF SDCRIRNGEN GIHFLNIRDG KRRGDALIEE ESEQDVQKAL EKHRMYMGQR YVEVYEINNE DVDALMKSLQ VKSSPVVNDG VVRLRGLPYS CNEKDIVDFF AGLNIVDITF VMDYRGRRTK GEAYVQFEED EMANQALLKH REEIGNRYIE IFPSRRNEVR THVGSYKGKK IASFPTAKYI TEPEMVFEED EVNEDIQPMY AFESKEKIEL PKEVPEKLPE AADFGTTSSL HFVHMRGLPF QANAQDIINF FAPLKPVRIT MEYSSSGKAT GEADVHFETH EDAVAAMLKD RSHVHHRYIE LFLNSCPKGG

image



Background

Alternative Names G-rich sequence factor 1 (GRSF-1)

Product Name: Recombinant Human GRSF1 Protein
Catalog #: PEH8528



Note

For Research Use Only , Not for Diagnostic Use