


Product Name: Recombinant Human CRYGC Protein
Catalog #: PEH8529



Summary

Name	CRGC_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P07315
Host	E.coli
Species	Human
Molecular Mass	22 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MGKITFYEDR AFQGRSYETT TDCPNLQPYF SRCNSIRVES GCWMLYERPN YQGQQYLLRR GEYDPYQQWM GLSDSIRSCC LIPQTVSHRL RLYEREDHKG LMMELSEDCP SIQDRFHLSE IRSLHVLEGC WVLYELPNYR GRQYLLRPQE YRRCQDWGAM DAKAGSLRRV VDLY

image



Background

Alternative Names	Gamma-crystallin C (Gamma-C-crystallin) (Gamma-crystallin 2-1) (Gamma-crystallin 3)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human CRYGC Protein
Catalog #: PEH8529

