


Product Name: Recombinant Human DEPDC5 Protein
Catalog #: PEH8531



Summary

Name	DEPD5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & N-His
Accession #	O75140
Host	E.coli
Species	Human
Molecular Mass	35 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	YPVLRLEQA EGFPQGDNST SAQGNYLEAI NLSFNVFDKH YINRNFDRGT QMSVVITPGV GVFEVDRLLM ILTKQRMIDN GIGVDLVCMD EQPLHAVPLF KLHNRSAPRD SRLGDDYNIP HWINHSFYTS KSQLFCNSFT PRIKLAGKKP ASEKAKNGRD TSLGSPKESE NALPIQVDYD

image



Background

Alternative Names GATOR complex protein DEPDC5 (DEP domain-containing protein 5)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human DEPDC5 Protein
Catalog #: PEH8531

