


Product Name: Recombinant Human NPRL3 Protein
Catalog #: PEH8532



Summary

Name	NPRL3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q12980
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PRPREDDVPF TARVGGRSL TPNALSFGSP TSSDDMTLTS PSMDNSSAEL LPSGDSPLNQ RMTENLLASL SEHERAAILS VPAAQNPEDL

image



Background

Alternative Names	GATOR complex protein NPRL3 (-14 gene protein) (Alpha-globin regulatory element-containing gene protein) (Nitrogen permease regulator 3-like protein) (Protein CGTHBA)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use