



Product Name: Recombinant Human GLCCI1 Protein
Catalog #: PEH8534



Summary

Name	GLCI1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q86VQ1
Host	E.coli
Species	Human
Molecular Mass	40 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PPAAAAPAEQ APRAKGRPRR SPESHRRSSS PERRSPGSPV CRADKAKSQQ VRTSSTIRRT SSLDTITGPY LTGQWPRDPH VHYPSCMKDK ATQTPSCWAE EGAEKRSHQR SASWGSADQL KEQIAKLRQQ LQRSKQSSRH SKEKDRQSPL HGNHITISHT QATGSRSVPM PLSNISVPKS SVSRVPCNVE GISPELEKVF IKENNGKEEV SKPLDIPDGR RAPLPAHYRS SSTRSIDTQT PSVQERSSSC SSHSPCVSPF CPPEQDQSP CSTEDLLYDR DKDSGSSSPL PKYASSPKPN NSYMFKREPP EGCERVKVFE EMASRQPISA PLFSCPDKNK

image



Background

Alternative Names Glucocorticoid-induced transcript 1 protein

Product Name: Recombinant Human GLCCI1 Protein
Catalog #: PEH8534



Note

For Research Use Only , Not for Diagnostic Use