


Product Name: Recombinant Human GOLGA4 Protein
Catalog #: PEH8537



Summary

Name	GOGA4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q13439
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MEAQHNDLEF KLAGAEREKQ KLGKEIVRLQ KDLRMLRKEH QQELEILKKE YDQEREKIK QEQEDLELKH NSTLKQLMRE FNTQLAQKEQ ELEMENTIKETI NKAQEVEAEL LESHQEETNQ LLKKIAEKDD DLKRTAKRYE EILDAREEEM TAKVRDLQTQ

image



Background

Alternative Names Golgin subfamily A member 4 (256 kDa golgin) (Golgin-245) (Protein 72.1) (Trans-Golgi p230)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human GOLGA4 Protein
Catalog #: PEH8537

