


Summary

Name	GGYF2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6Y7W6
Host	E.coli
Species	Human
Molecular Mass	45 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PTVWEGGSVW DLPLDTTTPG PALEQLQQL E KAKAAKLEQE RREAEMRAKR EEEEERKRQEE LRRQQEEILR RQQEERKRR EEEELARRKQ EEALRRQREQ EIALRRQREE EERQQQEEAL RRLEERRREE EERRKQEELL RKQEEEEAAKW AREEEEAQRR LEENRLRMEE EAARLRHEEE ERKRKELEVQ RQKELMRQRQ QQQEALRRLQ QQQQQQLAQ MKLPSSTWG QQSNTTACQS QATLSLAEIQ KLEERERQL REEQRRQRE LMKALQQQQQ QQQQKLSGWG NVSKPSGTTK SLLEIQQEEA RQMQRQQQQQ QQHQQPNRAR NNTHSNLHTS IGNSVWGSIN

image



Background

Alternative Names	GRB10-interacting GYF protein 2 (PERQ amino acid-rich with GYF domain-containing protein 2) (Trinucleotide repeat-containing gene 15 protein)
--------------------------	---

Product Name: Recombinant Human GIGYF2 Protein
Catalog #: PEH8538



Note

For Research Use Only , Not for Diagnostic Use