


**Product Name: Recombinant Human GNAT3 Protein**  
**Catalog #: PEH8539**



## Summary

<b>Name</b>	GNAT3_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	A8MTJ3
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	41 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 12 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	ESAKRSKELE KKLQEDAERD ARTVKLLLLG AGESGKSTIV KQMKIIHKNG YSEQECMEFK AVIYSNTLQS ILAIVKAMTT LGIDYVNPRS AEDQRQLYAM ANTLEDGGMT PQLAEVIKRL WRDPGIQACF ERASEYQLND SAAYYLNDLD RITASGYVPN EQDVLHSRVK TTGIIETQFS FKDLHFRMFD VGGQRSERKK WIHCFEGVTC IIFCAALSAY DMVLVEDEEV NRMHESLHLF NSICNHKYFS TTSIVLFLNK KDIFQEKVTK VHLSICFPEY TGPNTFEDAG NYIKNQFLDL NLKKEDKEIY SHMTCATDTQ NVKFVFDVAVT DIIKENLKD

## image



## Background

<b>Alternative Names</b>	Guanine nucleotide-binding protein G(t) subunit alpha-3 (Gustducin alpha-3 chain)
--------------------------	---

**Product Name: Recombinant Human GNAT3 Protein**  
**Catalog #: PEH8539**

---



**Note**

For Research Use Only , Not for Diagnostic Use