


Product Name: Recombinant Human HSPA14 Protein
Catalog #: PEH8542



Summary

Name	HSP7E_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q0VDF9
Host	E.coli
Species	Human
Molecular Mass	40 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EIDTGEETKF VNPEDVARLI FSKMKETAHS VLGSDANDVV ITVPDFGEEK QKNALGEAAR AAGFNVRLRI HEP SAALLAY GIGQDSPTGK SNILVFKLGG TSLSLVMEV NSGIYRVLST NTDDNIGGAH FTETLAQYLA SEFQRSFKHD VRGNARAMMK LTNSAEVAKH SLSTLGSANC FLDSLYEGQD FDCNVSRRARF ELLCSPFNK CIEAIRGLLD QNGFTADDIN KVVLCGGSSR IPKLQQLIKD LFP AVELLNS IPPDEVIPIG AAIEAGILIG KENLLVEDSL MIECSARDIL VKGVDESGAS RFTVLFPSGT PLPARRQHTL QAPGSISSVC LELYESDGKN

image



Background

Alternative Names Heat shock 70 kDa protein 14 (HSP70-like protein 1) (Heat shock protein HSP60)

Product Name: Recombinant Human HSPA14 Protein
Catalog #: PEH8542



Note

For Research Use Only , Not for Diagnostic Use