


Product Name: Recombinant Human IQSEC2 Protein
Catalog #: PEH8545



Summary

| | |
|------------------------------|---|
| Name | IQEC2_HUMAN |
| Purity | Greater than 85% as determined by reducing SDS-PAGE |
| Endotoxin level | Please contact with the lab for this information |
| Tag | C-His |
| Accession # | Q5JU85 |
| Host | E.coli |
| Species | Human |
| Molecular Mass | 33 kDa |
| Buffer | Phosphate buffered saline |
| Shipping | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below. |
| Stability&Storage | Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| Sequence | HSTSTSTSPA TTLQRKSDGE NSRTVSVEGD APGSDLSTAV DSPGSQPPYR LSQLPPSSSH MGGPPAGVGL PWAQRARLQP ASVALRKQEE EEIKRSKALS DSYELSTD LQ DKKVEMLERK YGGSFLSRRRA ARTIQTAFRQ YRMNKNFERL RSSASESRMS RRIILSNMRM QFSFEEYEKA QNPAYFEGKP ASLDEGAMAG ARSHRLERGL PYGGSCGGGI DGGGSSVTTS GEFSNDITEL EDSFSKQVKS LAESIDEALN CHPSGPMSEE PGSAQLEKRE SKEQQEDSSA |

image



Background

Alternative Names IQ motif and SEC7 domain-containing protein 2

Product Name: Recombinant Human IQSEC2 Protein
Catalog #: PEH8545



Note

For Research Use Only , Not for Diagnostic Use