


Product Name: Recombinant Human KANK1 Protein
Catalog #: PEH8548



Summary

Name	KANK1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q14678
Host	E.coli
Species	Human
Molecular Mass	40 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TSEGKPISSL DAFPTQEGTL SPVNLTDQI AAGLYACTNN ESTLKSIMKK KDGNKDSNGA KKNLQFVGIN GGYETTSSDD SSSDESSSSE SDDECDVIEY PLEEEEEED EDTRGMAEGH HAVNIEGLKS ARVEDEMQVQ ECEPEKVEIR ERYELSEKML SACNLLKNTI NDPKALTSKD MRFLNLTQHQ EWFRVSSQKS AIPAMVGDYI AAFEAISPDV LRYVINLADG NGNTALHYSV SHSNFEIVKL LLDADVCNVD HQNKAGYTPI MLAALAAVEA EKDMRIVEEL FGCGDVNAKA SQAGQTALML AVSHGRIDMV KGLLACGADV NIQDDEGSTA LMCASEHGHV

image



Background

Alternative Names KN motif and ankyrin repeat domain-containing protein 1 (Ankyrin repeat domain-containing protein 15) (Kidney ankyrin repeat-containing protein)

Product Name: Recombinant Human KANK1 Protein
Catalog #: PEH8548



Note

For Research Use Only , Not for Diagnostic Use