

Product Name: Recombinant Human LRRK1 Protein
Catalog #: PEH8550



Summary

Name	LRRK1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q38SD2
Host	E.coli
Species	Human
Molecular Mass	53 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SHLRLPWVDL DWLIDISCQI TELDLSANCL ATLPSVIPWG LINLRKLNLS DNHLGELPGV QSSDEIICSR LLEIDISSNK LSHLPPGFLH LSKLQKLTAS KNCLEKLFEE ENATNWIGLR KLQELDISDN KLTELPALFL HSFKSLNSLN VSRNNLKVFP DPWACPLKCC KASRNAECL PDKMAVFWKN HLKDVDLSEN ALKEVPLGLF QLDALMFLRL QGNQLAALPP QEKWTCRQLK TLDLSRNQLG KNEDGLKTKR IAFFTTRGRQ RSGTEAASVL EFPAPLSESL EVLCLNDNHL DTVPPSVCLL KSLSELYLGN NPGLRELPE LGQLGNLWQL DTEDLTISNV

image



SDS



WB

Background

Alternative Names Leucine-rich repeat serine/threonine-protein kinase 1 (EC 2.7.11.1)

Product Name: Recombinant Human LRRK1 Protein
Catalog #: PEH8550



Note

For Research Use Only , Not for Diagnostic Use