


Product Name: Recombinant Human TIMM44 Protein
Catalog #: PEH8551



Summary

Name	TIM44_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O43615
Host	E.coli
Species	Human
Molecular Mass	28 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LEESDVLQEA RRKYKTIESE TVRTSEVLRK KLGELTGTVK ESLHEVSKSD LGRKIKEGVE EAAKTAKQSA ESVSKGGEKL GRTAAFRALS QGVESVKKEI DDSVLGQTGP YRRPQRLRKR TEFAGDKFKE EKVFEPNEEA LGVVLHKDSK WYQQWKDFKE NNVVFNRFEE MKMKYDESDN AFIRASRALT DKVTDLLGGL FSKTEMSEVL TEILRVDPAF DKDRFLKQCE

image



Background

Alternative Names Mitochondrial import inner membrane translocase subunit TIM44

Note

Product Name: Recombinant Human TIMM44 Protein
Catalog #: PEH8551



For Research Use Only , Not for Diagnostic Use