


Product Name: Recombinant Human ZFP36L2 Protein
Catalog #: PEH8553



Summary

Name	TISD_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P47974
Host	E.coli
Species	Human
Molecular Mass	13 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ALLNKENKFR DRSFSENGDR SQHLLHLQQQ QKGGGGSQIN STRYKTELCR PFEESGTCKY GEKCQFAHGF HELRSLTRHP KYKTELCRTF

image



Background

Alternative Names	mRNA decay activator protein ZFP36L2 (Butyrate response factor 2) (EGF-response factor 2) (ERF-2) (TPA-induced sequence 11d) (Zinc finger protein 36, C3H1 type-like 2) (ZFP36-like 2)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use