


Product Name: Recombinant Human NOXA1 Protein
Catalog #: PEH8558



Summary

| | |
|------------------------------|---|
| Name | NOXA1_HUMAN |
| Purity | Greater than 85% as determined by reducing SDS-PAGE |
| Endotoxin level | Please contact with the lab for this information |
| Tag | C-His |
| Accession # | Q86UR1 |
| Host | E.coli |
| Species | Human |
| Molecular Mass | 10 kDa |
| Buffer | Phosphate buffered saline |
| Shipping | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below. |
| Stability&Storage | Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| Sequence | GKAKVVASAI PDDQGWGVRP QQPQGPGANH DARSLIMDSP RAGTHQGPLD AETEVGADRC TSTAYQEQRQ QVEQVGKQAP |

image



Background

Alternative Names NADPH oxidase activator 1 (NOX activator 1) (Antigen NY-CO-31) (NCF2-like protein) (P67phox-like factor) (p51-nox)

Note

For Research Use Only , Not for Diagnostic Use