


Product Name: Recombinant Human NSMF Protein
Catalog #: PEH8559



Summary

Name	NSMF_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6X4W1
Host	E.coli
Species	Human
Molecular Mass	25 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LQSSHCTLDE AFEDLDWDTE KGLEAVACDT EGFVPPKVML ISSKVPKAEY IPTIIRDDP SIIPILYDHE HATFEDILEE IERKLVYHK GAKIWKMLIF CQGGPGHLYL LKNKVATFAK VEKEEDMIHF WKRLSRLMSK VNPEPNVIHI MGCYILGNPN GEKLFQNLRT LMTPYRVTFE SPLELSAQGK QMIETYFDFR

image



Background

Alternative Names	NMDA receptor synaptonuclear signaling and neuronal migration factor (Nasal embryonic luteinizing hormone-releasing hormone factor) (Nasal embryonic LHRH factor)
--------------------------	---

Note

Product Name: Recombinant Human NSMF Protein
Catalog #: PEH8559



For Research Use Only , Not for Diagnostic Use