



Product Name: Recombinant Human NUP107 Protein
Catalog #: PEH8562



Summary

Name	NU107_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P57740
Host	E.coli
Species	Human
Molecular Mass	17 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SVTEFEQRHH CLELAKEADL DVATITKTVV ENIRKKDNGE FSHHDLAPAL DTGTTEEDRL KIDVIDWLVF DPAQRAEALK QGNAIMRKFL ASKKHEAAKE VFVKIPQDSI AEIYNQCEEQ GMESPLPAED

image



Background

Alternative Names Nuclear pore complex protein Nup107 (107 kDa nucleoporin) (Nucleoporin Nup107)

Note

For Research Use Only , Not for Diagnostic Use