

Product Name: Recombinant Human GLE1 Protein
Catalog #: PEH8567



Summary

Name	GLE1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q53GS7
Host	E.coli
Species	Human
Molecular Mass	31 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MKGTEGLRLW QEEQERKVQA LSEMASEQLK RFDEWKELKQ HKEFQDLREV MEKSSREALG HQEKLKAEHR HRAKILNLKL REAEQQRVKQ AEQERLRKEE GQIRLRALYA LQEEMLQLSQ QLDASEQHKA LLKVDLAAFQ TRGNQLCSLI SGIIRASSES SYPTAESQAE AERALREMRD LLMNLGQEIT RACEDKRRQD EEEAQVKLQE AQMQQGPEAH KEPPAPSQGP GKGQNEDLQV KVQDITMQWY

image



Background

Alternative Names Nucleoporin GLE1 (hGLE1) (GLE1-like protein)

Note

Product Name: Recombinant Human GLE1 Protein
Catalog #: PEH8567



For Research Use Only , Not for Diagnostic Use