


Product Name: Recombinant Human EHHADH Protein
Catalog #: PEH8568



Summary

Name	ECHP_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q08426
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	NQLATANKMI TSVLEKEASK MQQSGHPWSG PKPRLTSSVK ELGGVDLVIE AVFEEMSLKK QVFAELSAVC KPEAFLCTNT SALDVDEIAS

image



Background

Alternative Names	Peroxisomal bifunctional enzyme (PBE) (PBFE) (L-bifunctional protein) (LBP) (Multifunctional enzyme 1) (MFE1) [Includes: Enoyl-CoA hydratase/3,2-trans-enoyl-CoA isomerase (EC 4.2.1.17) (EC 5.3.3.8); 3-hydroxyacyl-CoA dehydrogenase (EC 1.1.1.35)]
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human EHHADH Protein
Catalog #: PEH8568

