


Product Name: Recombinant Human PEX19 Protein
Catalog #: PEH8569



Summary

Name	PEX19_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P40855
Host	E.coli
Species	Human
Molecular Mass	12 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	FQELFDSELA SQATAEFEKA MKELAESEPH LVEQFQKLSE AAGRVGSDMT SQQEFTSCLK ETLGSLAKNA TDLQNSSMSE EELTKAMEGL

image



Background

Alternative Names	Peroxisomal biogenesis factor 19 (33 kDa housekeeping protein) (Peroxin-19) (Peroxisomal farnesylated protein)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use