

Product Name: Recombinant Human BFSP2 Protein
Catalog #: PEH8570



Summary

Name	BFSP2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q13515
Host	E.coli
Species	Human
Molecular Mass	37 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LGGCLVEYMA KVHALEQVSQ ELETQLRMHL ESKATRSQNW GALRASWASS CQQVGEAVLE NARLMLQTET IQAGADDFKE RYENEQPFK AEEEEINSLY KVIDEANLTK MDLESQIESL KEELGSLSRN YEEDVKLLHK QLAGCELEQM DAPIGTGLDD ILETIRIQWE RDVEKNRVEA GALLOAKQQA EVAHMSQTQE EKLAALRVE LHNTSCQVQS LQAETESLRA LKRGLENTLH DAKHWHDMEL QNLGAVVGRL EAELEIRAE AEQQQQRH LLARKCQLQK DVASYHALLD

image



Background

Alternative Names	Phakinin (49 kDa cytoskeletal protein) (Beaded filament structural protein 2) (Lens fiber cell beaded filament protein CP 47) (CP47) (Lens fiber cell beaded filament protein CP 49) (CP49) (Lens intermediate filament-like light) (LIFL-L)
--------------------------	--

Product Name: Recombinant Human BFSP2 Protein
Catalog #: PEH8570



Note

For Research Use Only , Not for Diagnostic Use