


Product Name: Recombinant Human PHLDB1 Protein
Catalog #: PEH8572



Summary

Name	PHLB1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q86UU1
Host	E.coli
Species	Human
Molecular Mass	40 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VTRERKNSIT EISDNEDLL EYHRRQRQER LREQEMERLE RQRLETILNL CAEYSRADGG PEAGELPSIG EATAALALAG RRPSRGLAGA SGRSSEEPGV ATQRLWESME RSDEENLKEE CSSTESTQQE HEDAPSTKLQ GEVLALEEEER AQVLGHVEQL KVRVKELEQQ LQESAREAEM ERALLQGERE AERALLQKEQ KAVDQLQEKL VALETGIQKE RDKEAEALET ETKLFEDLEF QQLERESRVE EERELAGQGL LRSKAELLRS IAKRKERLAI LDSQAGQIRA QAVQESERLA RDKNASLQLL QKEKEKLTVL ERRYHSLTGG

image



Background

Alternative Names Pleckstrin homology-like domain family B member 1 (Protein LL5-alpha)

Product Name: Recombinant Human PHLDB1 Protein
Catalog #: PEH8572



Note

For Research Use Only , Not for Diagnostic Use