


Product Name: Recombinant Human PCBP4 Protein
Catalog #: PEH8573



Summary

Name	PCBP4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P57723
Host	E.coli
Species	Human
Molecular Mass	31 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EPELSITLTL RMLMHGKEVG SIIGKKGETV KRIREQSSAR ITISEGSCPE RITTITGSTA AVFHAVSMIA FKLDEDLCAA PANGGNVSRP PVTLRVIPA SQCGSLIGKA GTKIKEIRET TGAQVQVAGD LLPNSTERAV TVSGVPDAII LCVRQICAVI LESPPKGATI PYHPSLSLGT VLLSANQGFS VQGQYGAVTP AEVTKLQQLS SHAVPFATPS VVPGLDPGTQ TSSQEFLVPN DLIGCVIGRQ GSKISEIRQM SGAHKIGNQ AEGAGERHVT

image



Background

Alternative Names Poly(rC)-binding protein 4 (Alpha-CP4)

Product Name: Recombinant Human PCBP4 Protein
Catalog #: PEH8573



Note

For Research Use Only , Not for Diagnostic Use