


Product Name: Recombinant Human POU4F3 Protein
Catalog #: PEH8574



Summary

Name	PO4F3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q15319
Host	E.coli
Species	Human
Molecular Mass	12 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ESDPRELEAF AERFKQRRIK LGVTQADVGA ALANLKIPGV GSLSQSTICR FESLTLSHNN MIALKPVLQA WLEEAEAAAYR EKNSKPELFN

image



Background

Alternative Names	POU domain, class 4, transcription factor 3 (Brain-specific homeobox/POU domain protein 3C) (Brain-3C) (Brn-3C)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use