

Product Name: Recombinant Human DTX2 Protein
Catalog #: PEH8575



Summary

Name	DTX2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q86UW9
Host	E.coli
Species	Human
Molecular Mass	29 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	HGSVKRLRKM SVKGATPKPE PEPEQVIKNY TEELKVPPDE DCIICMEKLS TASGYSDVTD SKAIGSLAVG HLTKCSHAFH LLCLLAMYCN GNKDGSLQCP SCKTIYGEKT GTQPQ GKMEV LRFQMSLP GH EDCGTILIVY SIPHGIQGPE HPNPGKPFTA RGFPRQCYLP DNAQGRKVLE LLKVAWKRRLL IFTVGTSSST GETDVTVVWNE IHHKTEMDRN ITGHGYDPDN YLQNVLAELA

image



Background

Alternative Names	Probable E3 ubiquitin-protein ligase DTX2 (EC 2.3.2.27) (Protein deltex-2) (Deltex2) (hDTX2) (RING finger protein 58) (RING-type E3 ubiquitin transferase DTX2)
--------------------------	---

Product Name: Recombinant Human DTX2 Protein
Catalog #: PEH8575



Note

For Research Use Only , Not for Diagnostic Use