


Product Name: Recombinant Human PRMT2 Protein
Catalog #: PEH8580



Summary

Name	ANM2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P55345
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ESQGEEPAEC SEAGLLQEGV QPEEFVAIAD YAATDETQLS FLRGEKILIL RQTTADWWWG ERAGCCGYIP

image



Background

Alternative Names	Protein arginine N-methyltransferase 2 (EC 2.1.1.319) (Histone-arginine N-methyltransferase PRMT2)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use