

Product Name: Recombinant Human FAM13A Protein
Catalog #: PEH8581




Summary

Name	FA13A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O94988
Host	E.coli
Species	Human
Molecular Mass	43 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SSTKLSLHD NQDGLVNMES LNSTRSHERT GPDDFEWMSD ERKGNEKDGG HTQHFESPTM KIQEHPSLSD TKQQRNQDAG DQEESFVSEV PQSDLTALCD EKNWEEPIPA FSSWQRENSD SDEAHLSPQA GRLIRQLLDE DSDPMLSPRF YAYGQSRQYL DDTEVPPSP NSHSFMRRRS SSLGSYDDEQ EDLTPAQLTR RIQSLKKKIR KFEDRFEEK KYRPSHSDKA ANPEVLKWTN DLAKFRRQLK ESKLKISEED LTPRMRQSN TLPKSFQSGL EKEDEKKQEL VDKAIKPSVE ATLESIQRKL QEKRAESSRP EDIKDMTKDQ IANEKVALQK ALLYYESIHG

image



SDS



WB

Background

Alternative Names Protein FAM13A

Product Name: Recombinant Human FAM13A Protein
Catalog #: PEH8581



Note

For Research Use Only , Not for Diagnostic Use