


Product Name: Recombinant Human FAM161A Protein
Catalog #: PEH8582



Summary

Name	F161A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q3B820
Host	E.coli
Species	Human
Molecular Mass	31 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	HLQNSSLPC RSACGCRNPR CPEQAVKLKC KHKVRCPTPD FEDLPERYQK HLSEHKSPKL LTVCKPFDLH ASPHASIKRE KILADIEADE ENLKETRWPY LSPRRKSPVR CAGVNPVPCN CNPPVPTVSS RGREQAVRKS EKERMREYQR ELEEREKLLK KRPLLFERVA QKNARMAAEK HYSNTLKALG ISDEFVSKKG QSGKVLEYFN NQETKSVTED KESFNEEEEI EERENGEEY FIDTNSQDSY KEKDEANES

image



Background

Alternative Names Protein FAM161A

Product Name: Recombinant Human FAM161A Protein
Catalog #: PEH8582



Note

For Research Use Only , Not for Diagnostic Use