


Product Name: Recombinant Human MTSS2 Protein
Catalog #: PEH8586



Summary

Name	MTSS2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q765P7
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ASETCQSVSE CSSPTSDWSK VGSHEQPSGA TLQRRKDRVE LLRDTEPGPA SGGTLGPSGE EAPRPRMSPA TIAAKHGEEV SPAASDLAMV LTRGLSLEHQ KSSRDSLQYS SGYSTQTTTP SCSED TIPSQ GSDYDCYSVN GDADEGPPE FDKSSTIPRN SNIAQNYRRL IQTKRPASTA GLPTAGLPTA TGLPSGAPPG

image



Background

Alternative Names Protein MTSS 2 (Actin-bundling with BAIAP2 homology protein 1) (ABBA-1) (MTSS1-like protein)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human MTSS2 Protein
Catalog #: PEH8586

