


Product Name: Recombinant Human PID1 Protein
Catalog #: PEH8589



Summary

Name	PCL11_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q7Z2X4
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	CKVTYLGVKS TTGMQFLSGC TEKPVIELWK KHTLAREDVF PANALLEIRP FQVWLHHLDH KGEATVHMDT FQVARIAYCT ADHNVSPNIF AWVYREINDD LSYQMDCHAV ECESKLEAKK LAHAMMEAFR KTFHSMKSDG RIHSNSSEE VSQELESDDG

image



Background

Alternative Names	PTB-containing, cubilin and LRP1-interacting protein (P-CL11) (Phosphotyrosine interaction domain-containing protein 1) (Protein NYGGF4)
--------------------------	--

Note

Product Name: Recombinant Human PID1 Protein
Catalog #: PEH8589



For Research Use Only , Not for Diagnostic Use