

**Product Name: Recombinant Human PLEKHG4 Protein**  
**Catalog #: PEH8590**



## Summary

<b>Name</b>	PKHG4_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	Q58EX7
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	20 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	SLAIKQAWTA DISHLLWRQA VHNKEVRMAE MVSMGVGNKA FRDIAPSEEA INDRTVNYVL KCREVRSRAS IAVAPFDHDS LYLGASNSLP GDPASCSVLG SLNLHLYRDP ALLGLRCPLY PSFP EEAALE AEAELGGQPS LTAEDSEISS QCPSASGSSG SDSSCVSGQA

## image



## Background

<b>Alternative Names</b>	Puratrophin-1 (Pleckstrin homology domain-containing family G member 4) (PH domain-containing family G member 4) (Purkinje cell atrophy-associated protein 1)
--------------------------	---

## Note

**Product Name: Recombinant Human PLEKHG4 Protein**  
**Catalog #: PEH8590**



---

For Research Use Only , Not for Diagnostic Use