


**Product Name: Recombinant Human RAPGEF3 Protein**  
**Catalog #: PEH8593**



## Summary

<b>Name</b>	RPGF3_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	O95398
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	18 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	LAEAVALLSQ RGPDALLTVA LRKPPGQRTD EELDIFEEL LHIKAVAHLS NSVKRELAHV LLFEPHSKAG TVLFSQGDKG TSWYIIWKGS VNVVTHGKGL VTTLHEGDDF GQLALVNDAP RAATIILRED NCHFLRVDKQ DFNRIIKDVE

## image



## Background

<b>Alternative Names</b>	Rap guanine nucleotide exchange factor 3 (Exchange factor directly activated by cAMP 1) (Exchange protein directly activated by cAMP 1) (EPAC 1) (Rap1 guanine-nucleotide-exchange factor directly activated by cAMP) (cAMP-regulated guanine nucleotide exchange factor I) (cAMP-GEFI)
--------------------------	---

## Note

**Product Name: Recombinant Human RAPGEF3 Protein**  
**Catalog #: PEH8593**



---

For Research Use Only , Not for Diagnostic Use