


Product Name: Recombinant Human RIN1 Protein
Catalog #: PEH8594



Summary

Name	RIN1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q13671
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	THCFQHLLRV AYQDPSSGCT SKTLAVPPEA SIATLNQLCA TKFRVTQPNT FGLFLYKEQG YHRLPPGALA HRLPTTGYLV YRRAEWPETQ GAVTEEEGSG QSEARSRGEE QGCQGDGDAG VKASPRDIRE QSETTAEGGQ GQAQEGPAQP GEPEAECSRA

image



Background

Alternative Names	Ras and Rab interactor 1 (Ras inhibitor JC99) (Ras interaction/interference protein 1)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human RIN1 Protein
Catalog #: PEH8594

