


Product Name: Recombinant Human RASGRF2 Protein
Catalog #: PEH8599



Summary

Name	RGRF2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O14827
Host	E.coli
Species	Human
Molecular Mass	17 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QSPAASPPPH TGQIPLDSR GLSSPEQSPG TVEENVDPNR VDLCKNLKRS IQKAVLESAP ADRAVVESSP AADTTELSPC RSPSTPRHLR YRQPGGQTAD NAHCSVSPAS AFAIATAAAG HGSPPGFNNT ERTCDKEFII

image



Background

Alternative Names	Ras-specific guanine nucleotide-releasing factor 2 (Ras-GRF2) (Ras guanine nucleotide exchange factor 2)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use