


Summary

Name	RGS12_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O14924
Host	E.coli
Species	Human
Molecular Mass	41 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RSVEVARGRA GYGFTLSGQA PCVLSCVMRG SPADFVGLRA GDQILAVNEI NVKKASHEDV VKLIGKCSGV LHMVIAEGVG RFESCSSDEE GGLYEGKGWL KPKLDSKALG INRAERVVEE MQSGGIFNMI FENPSLCASN SEPLKCLKQRS LSESAATRFD VGHESINNPV PNMLSKEEIS KVIHDDSVFS IGLSHDDFA LDASILNVAM IVGYLGSIEL PSTSSNLESD SLQAIRGCMR RLRAEQKIHS LVTMIMHDC VQLSTDKAGV VAEYPAEKLA FSAVCPDDRR FFGLVMTMQTN DDGSLAQEEE GALRTSCHVF MVDPDFNHNK IHQGIARRFG FECTADPDTN GCLEFPASSL

image



Background

Alternative Names Regulator of G-protein signaling 12 (RGS12)

Product Name: Recombinant Human RGS12 Protein
Catalog #: PEH8600



Note

For Research Use Only , Not for Diagnostic Use