


Product Name: Recombinant Human ARHGAP21 Protein
Catalog #: PEH8603



Summary

Name	RHG21_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q5T5U3
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PVNHLLHQNSL LNQQTWVRTD SAPDQQVETG KPSLSGASA KPAPQSSENA GTSDELPLVS QRNQDLSLQE AETEQSDTLD NKEAVILREK PPSGRQTPQP LRHQSYILAV NDQETGSDTT CWLPNDARRE VHIKRMEEK ASSTSPPGDS LASIPFIDEP TSPSIDHDIA HIPASAVISA STSQVPSIAT

image



Background

Alternative Names	Rho GTPase-activating protein 21 (Rho GTPase-activating protein 10) (Rho-type GTPase-activating protein 21)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human ARHGAP21 Protein
Catalog #: PEH8603

