


Product Name: Recombinant Human ARHGAP4 Protein
Catalog #: PEH8604



Summary

Name	RHG04_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P98171
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DGDEVAEICV EMELRDEILP RAQNIQSRLD RQTETEEVN KTLKATLQAL LEVVASDDGD VLDSFQTSPS TESLKSTSSD

image



Background

Alternative Names	Rho GTPase-activating protein 4 (Rho-GAP hematopoietic protein C1) (Rho-type GTPase-activating protein 4) (p115)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use