

Product Name: Recombinant Human ARHGEF6 Protein
Catalog #: PEH8606



Summary

Name	ARHG6_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q15052
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EPPQIIKPWS LSCLRPAPPL RPSAALGYKE RMSYILKESS KSPKTMKKFL HKRKTERKPS EEEYVIRKST AALEEDAQIL KVIEAYCTSA NFQQGHGSST RKDSIPQVLL PEEKLIIEE TRSNGQTIME EKSLVDTVYA LKDEVRELKQ ENKRMKQCLE

image



Background

Alternative Names	Rho guanine nucleotide exchange factor 6 (Alpha-Pix) (COOL-2) (PAK-interacting exchange factor alpha) (Rac/Cdc42 guanine nucleotide exchange factor 6)
--------------------------	--

Note

Product Name: Recombinant Human ARHGEF6 Protein
Catalog #: PEH8606



For Research Use Only , Not for Diagnostic Use