


Summary

Name	SPB7_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	O75635
Host	E.coli
Species	Human
Molecular Mass	76 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SSLSLFAALA LURLGAQDDS LSQIDKLLHV NTASGYGNSS NSQSGLQSQL KRVFSDINAS HKDYDLSIVN GLFAEKVYGF HKDYIECAEK LYDAKVERVD FTNHLEDTRR NINKWVENET HGKIKNVIGE GGISSAVMV LVNAVYFKGK WQSAFTKSET INCHFSPKC SGKAVAMMHQ ERKFNLVIE DPSMKILELR YNGGINMYVL LPENDLSEIE NKLTFQNLME WTNPRRMTSK YVEVFFPQFK IEKNYEMKQY LRALGLKDIF DESKADLSGI

image



Background

Alternative Names Serpin B7 (Megsin) (TP55)

Product Name: Recombinant Human SERPINB7 Protein
Catalog #: PEH8613



Note

For Research Use Only , Not for Diagnostic Use