


## Summary

<b>Name</b>	SPB7_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	N-MBP & C-His
<b>Accession #</b>	O75635
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	80 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	DNQGNGNVFF SLSLFAALA LVRLGAQDDS LSQIDKLLHV NTASGYGNSS NSQSGLQSQL KRVFSDINAS HKDYDLSIVN GLFAEKVYGF HKDYIECAEK LYDAKVERVD FTNHLEDTRR NINKWVENET HGKIKNVIGE GGISSAVMV LVNAVYFKGK WQSAFTKSET INCHFSPKC SGKAVAMMHQ ERKFNLVIE DPSMKILELR YNGGINMYVL LPENDLSEIE NKLTFQNLME WTNPRRMTSK YVEVFFPQFK IEKNYEMKQY LRALGLKDIF DESKADLSGI ASGGRLYISR MMHKSYIEVT EEGTEATAAT GSNIVEKQLP

## image



## Background

**Alternative Names** Serpin B7 (Megsin) (TP55)

**Product Name: Recombinant Human SERPINB7 Protein**  
**Catalog #: PEH8615**



---

**Note**

For Research Use Only , Not for Diagnostic Use