


Product Name: Recombinant Human STRA8 Protein
Catalog #: PEH8621



Summary

Name	STRA8_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q7Z7C7
Host	E.coli
Species	Human
Molecular Mass	28 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	NKAKSHIPEL EQTLDNLLKL KASFNLEDGH ASSLEEVKKE YASMYSGNDS FPQNGSSPWY LNFYKQTMDL LTGSGIITPQ EAALPIVSAA ISHLWQNLSE ERKASLRQAW AQKHRGPATL AEACREPACA EGSVKDSGVD SQGASCSLVS TPEILFEDA FDVASFLLDKS EVPSTSSSSS VLASCPENP EEKFQLYMQI INFFKGLSCA NTQVKQEASF PVDEEMIMLQ CTETFDDEDL

image



Background

Alternative Names Stimulated by retinoic acid gene 8 protein homolog

Note

Product Name: Recombinant Human STRA8 Protein
Catalog #: PEH8621



For Research Use Only , Not for Diagnostic Use