


Product Name: Recombinant Human TBC1D8 Protein
Catalog #: PEH8623



Summary

Name	TBCD8_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O95759
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	NMDQLIEFKA FVSCDIMYN GEMNEKIKLL YRLHIPPALT ENDRDSQSPL RNPLLSTSRP LVFGKPNGDA VDYQKQLKQM IKDLAKEKDK TEKELPKMSQ REFIQFCKTL YSMFHEDPEE NDLYQAIATV TLLLLQIGEV GQRGSSSGSC SQECGEELRA SAPSPEDSVF ADTGKTPQDS QAFPEAAERD

image



Background

Alternative Names	TBC1 domain family member 8 (AD 3) (Vascular Rab-GAP/TBC-containing protein)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human TBC1D8 Protein
Catalog #: PEH8623

