


Summary

Name	DHB3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P37058
Host	E.coli
Species	Human
Molecular Mass	12 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GKAYSFELAK RGLNVVLISR TLEKLEAIAT EIERTTGRSV KIIQADFTKD DIYEHIKEKL AGLEIGILVN NVGMLPNLLP SHFLNAPDEI

image



Background

Alternative Names	Testosterone 17-beta-dehydrogenase 3 (EC 1.1.1.64) (17-beta-hydroxysteroid dehydrogenase type 3) (17-beta-HSD 3) (Estradiol 17-beta-dehydrogenase 2) (EC 1.1.1.62) (Short chain dehydrogenase/reductase family 12C member 2) (Testicular 17-beta-hydroxysteroid dehydrogenase)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human HSD17B3 Protein
Catalog #: PEH8625

