


Product Name: Recombinant Human SOX30 Protein
Catalog #: PEH8627



Summary

Name	SOX30_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O94993
Host	E.coli
Species	Human
Molecular Mass	35 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RASRVVKLEG PGPALGYFRG DEKGKLEAEE VMRDSMQGGA GKSPAAIREG VIKTEEPERL LEDCRLGAEP ASNGLVHGSA EVILAPTSGA FGPHQQDLRI PLTLHTVPPG ARIQFQGAPP SELIRLTKVP LTPVPTKMQS LLEPSVKIET KDVPLTVLPS DAGIPDTPFS KDRNGHVKRP MNAFMVWARI HRPALAKANP AANNAEISVQ LGLEWNKLSE EQKKPYYDEA QKIKEKHREE FPGWVYQPRP GKRKRFP LSV SNVFSGTTQN IISTNPTTVY PYRSPTYSVV IPSLQNPITH

image



Background

Alternative Names Transcription factor SOX-30

Product Name: Recombinant Human SOX30 Protein
Catalog #: PEH8627



Note

For Research Use Only , Not for Diagnostic Use