



Product Name: Recombinant Human TRO Protein
Catalog #: PEH8629



Summary

Name	TROP_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q12816
Host	E.coli
Species	Human
Molecular Mass	22 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TTEVTNTQAS SVTAQPKKAN KMKRVTAKEA QGSQSPTGHE GGTIQLKSPL QVLKLPVISQ NIHAPIANES ASSQALITSI KPKKASKAKK AANKAIASAT EVSLAATATH TATTQGQITN ETASIHHTAA SIRTCKASKA RKTIKAVINT DTEHIEALNV TDAATRQIEA SVVAIRPKKS KGKKAASRGP NSVSEISEAP

image



Background

Alternative Names Trophinin (MAGE-D3 antigen)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human TRO Protein
Catalog #: PEH8629

